

RP-Vantage™

An Innovative Remote Presence Solution for Surgical Telementoring and Collaboration

InTouch Health provides a suite of Remote Presence solutions that allows physicians to be available on demand whenever clinical expertise is needed. RP-Vantage™ is an FDA-cleared product designed to enable surgical telementoring and remote collaboration in ORs and procedure rooms.

With RP-Vantage, distant proctoring or learning surgeons can avoid the costs and inconveniences of travel and be virtually present to train and collaborate with on-site surgeons. Enabling more frequent surgical collaboration may contribute to improved clinical outcomes and accelerated adoption of innovative medical devices and procedures.



Key Features:

- Two-way audio/video communications
- 9' articulating boom arm with P-T-Z camera
- Multiple video-in ports integrating additional sources
- Live telestration and collaboration tools
- Video recording and image capture
- Reliable connectivity through SureConnect®
- Wi-Fi and battery enabled
- Cart-based platform that is mobile and shippable

Applications and Use Cases:

- Remote training via peer-to-peer proctoring and case observation
- Preoperative planning
- New procedure training/device proctoring
- Unscheduled procedures
- Clinical study site training and monitoring
- Clinical and technical product support and in-servicing

